



NAAA Antibody, FITC conjugated

Product Code	CSB-PA22919C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02083
Immunogen	Recombinant Human N-acylethanolamine-hydrolyzing acid amidase protein (29-199AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acylethanolamine-hydrolyzing acid amidase (EC 3.5.1.-) (Acid ceramidase-like protein) (N-acylsphingosine amidohydrolase-like) (ASAH-like protein) [Cleaved into: N-acylethanolamine-hydrolyzing acid amidase subunit alpha; N-acylethanolamine-hydrolyzing acid amidase subunit beta], NAAA, AS AHL PLT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAAA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.