



SH3BP1 Antibody, Biotin conjugated

Product Code	CSB-PA23069D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3L3
Immunogen	Recombinant Human SH3 domain-binding protein 1 protein (161-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding protein 1, SH3BP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SH3BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.