



GPR98 Antibody, HRP conjugated

Product Code	CSB-PA23089B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXG9
Immunogen	Recombinant Human G-protein coupled receptor 98 protein (6155-6306AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 98 (Monogenic audiogenic seizure susceptibility protein 1 homolog) (Usher syndrome type-2C protein) (Very large G-protein coupled receptor 1), GPR98, KIAA0686 KIAA1943 MASS1 VLGR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPR98
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.