



# POU5F1B Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA23129C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q06416
<b>Immunogen</b>	Recombinant Human Putative POU domain, class 5, transcription factor 1B protein (10-138AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Putative POU domain, class 5, transcription factor 1B (Oct4-pg1) (Octamer-binding protein 3-like) (Octamer-binding transcription factor 3-like), POU5F1B, OCT4PG1 OTF3C OTF3P1 POU5F1P1 POU5FLC20 POU5FLC8
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	POU5F1B
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.