



POU5F1B Antibody, Biotin conjugated

Product Code	CSB-PA23129D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06416
Immunogen	Recombinant Human Putative POU domain, class 5, transcription factor 1B protein (10-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative POU domain, class 5, transcription factor 1B (Oct4-pg1) (Octamer-binding protein 3-like) (Octamer-binding transcription factor 3-like), POU5F1B, OCT4PG1 OTF3C OTF3P1 POU5F1P1 POU5FLC20 POU5FLC8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POU5F1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.