



IGHM Antibody, FITC conjugated

Product Code	CSB-PA360788LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01871
Immunogen	Recombinant Human Immunoglobulin heavy constant mu protein (45-178AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin heavy constant mu (Ig mu chain C region) (Ig mu chain C region BOT) (Ig mu chain C region GAL) (Ig mu chain C region OU), IGHM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGHM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.