



KIR3DS1 Antibody, HRP conjugated

Product Code	CSB-PA614817LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14943
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 3DS1 protein (198-335AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 3DS1 (MHC class I NK cell receptor) (Natural killer-associated transcript 10) (NKAT-10), KIR3DS1, NKAT10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KIR3DS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.