



KIR3DS1 Antibody, FITC conjugated

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA614817LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14943 |
| Immunogen | Recombinant Human Killer cell immunoglobulin-like receptor 3DS1 protein (198-335AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Killer cell immunoglobulin-like receptor 3DS1 (MHC class I NK cell receptor) (Natural killer-associated transcript 10) (NKAT-10), KIR3DS1, NKAT10 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | KIR3DS1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |