



SGCB Antibody, HRP conjugated

Product Code	CSB-PA618082LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16585
Immunogen	Recombinant Human Beta-sarcoglycan protein (124-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-sarcoglycan (Beta-SG) (43 kDa dystrophin-associated glycoprotein) (43DAG) (A3b), SGCB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.