



ITGA7 Antibody, Biotin conjugated

Product Code	CSB-PA618797LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13683
Immunogen	Recombinant Human Integrin alpha-7 protein (175-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-7 [Cleaved into: Integrin alpha-7 heavy chain; Integrin alpha-7 light chain; Integrin alpha-7 70 kDa form], ITGA7
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ITGA7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.