



PPP2R5E Antibody, Biotin conjugated

Product Code	CSB-PA619075LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16537
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit epsilon isoform protein (2-149AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit epsilon isoform (PP2A B subunit isoform B'-epsilon) (PP2A B subunit isoform B56-epsilon) (PP2A B subunit isoform PR61-epsilon) (PP2A B subunit isoform R5-epsilon), PPP2R5E
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PPP2R5E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.