



ZP4 Antibody, Biotin conjugated

Product Code	CSB-PA621518LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12836
Immunogen	Recombinant Human Zona pellucida sperm-binding protein 4 protein (368-457AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zona pellucida sperm-binding protein 4 (Zona pellucida glycoprotein 4) (Zp-4) (Zona pellucida protein B) [Cleaved into: Processed zona pellucida sperm-binding protein 4], ZP4, ZPB
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ZP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.