



SNTB1 Antibody, Biotin conjugated

Product Code	CSB-PA621699LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13884
Immunogen	Recombinant Human Beta-1-syntrophin protein (51-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1-syntrophin (59 kDa dystrophin-associated protein A1 basic component 1) (DAPA1B) (BSYN2) (Syntrophin-2) (Tax interaction protein 43) (TIP-43), SNTB1, SNT2B1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNTB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.