



SEC23B Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA621874LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q15437 |
| Immunogen | Recombinant Human Protein transport protein Sec23B protein (208-325AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein transport protein Sec23B (hSec23B) (SEC23-related protein B), SEC23B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SEC23B |