



SEC23B Antibody, Biotin conjugated

Product Code	CSB-PA621874LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15437
Immunogen	Recombinant Human Protein transport protein Sec23B protein (208-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein transport protein Sec23B (hSec23B) (SEC23-related protein B), SEC23B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC23B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.