



SEMA3F Antibody, HRP conjugated

Product Code	CSB-PA622647LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13275
Immunogen	Recombinant Human Semaphorin-3F protein (565-651AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-3F (Sema III/F) (Semaphorin IV) (Sema IV), SEMA3F
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SEMA3F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.