



E2F5 Antibody, HRP conjugated

Product Code	CSB-PA624021LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15329
Immunogen	Recombinant Human Transcription factor E2F5 protein (193-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor E2F5 (E2F-5), E2F5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	E2F5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.