



LARP7 Antibody, HRP conjugated

Product Code	CSB-PA672326LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0J3
Immunogen	Recombinant Human La-related protein 7 protein (409-471AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	La-related protein 7 (La ribonucleoprotein domain family member 7) (P-TEFb-interaction protein for 7SK stability) (PIP7S), LARP7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LARP7