



LARP7 Antibody, HRP conjugated

Product Code	CSB-PA672326LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0J3
Immunogen	Recombinant Human La-related protein 7 protein (409-471AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	La-related protein 7 (La ribonucleoprotein domain family member 7) (P-TEFb-interaction protein for 7SK stability) (PIP7S), LARP7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LARP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.