



ANO6 Antibody, HRP conjugated

Product Code	CSB-PA679761LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4KMQ2
Immunogen	Recombinant Human Anoctamin-6 protein (12-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anoctamin-6 (Small-conductance calcium-activated nonselective cation channel) (SCAN channel) (Transmembrane protein 16F), ANO6, TMEM16F
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ANO6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.