



# ANO6 Antibody, HRP conjugated

Product Code	CSB-PA679761LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4KMQ2
Immunogen	Recombinant Human Anoctamin-6 protein (12-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anoctamin-6 (Small-conductance calcium-activated nonselective cation channel) (SCAN channel) (Transmembrane protein 16F), ANO6, TMEM16F
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ANO6