



ompA Antibody, HRP conjugated

Product Code	CSB-PA706579LB01RMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q52657
Immunogen	Recombinant Rickettsia conorii Outer membrane protein A protein (1734-2021AA)
Raised In	Rabbit
Species Reactivity	Rickettsia conorii
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein A (190 kDa antigen) (Cell surface antigen) (rOmp A) (rOmpA) [Cleaved into: 120 kDa surface-exposed protein (120 kDa outer membrane protein OmpA); 32 kDa beta peptide], ompA
Immunogen Species	Rickettsia conorii (strain ATCC VR-613 / Malish 7)
Research Area	Others
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.