



# YOD1 Antibody, HRP conjugated

Product Code	CSB-PA713352LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VVQ6
Immunogen	Recombinant Human Ubiquitin thioesterase OTU1 protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin thioesterase OTU1 (EC 3.4.19.12) (DUBA-8) (HIV-1-induced protease 7) (HIN-7) (HsHIN7) (OTU domain-containing protein 2), YOD1, DUBA8 HIN7 OTUD2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	YOD1