



YOD1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA713352LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5VVQ6 |
| Immunogen | Recombinant Human Ubiquitin thioesterase OTU1 protein (1-44AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ubiquitin thioesterase OTU1 (EC 3.4.19.12) (DUBA-8) (HIV-1-induced protease 7) (HIN-7) (HshIN7) (OTU domain-containing protein 2), YOD1, DUBA8 HIN7 OTUD2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | YOD1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |