



MAP3K21 Antibody, HRP conjugated

Product Code	CSB-PA719421LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCX8
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 21 protein (640-794AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 21 (EC 2.7.11.25) (Mitogen-activated protein kinase kinase kinase MLK4) (Mixed lineage kinase 4), MAP3K21, KIAA1804 MLK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K21