



MAP3K21 Antibody, FITC conjugated

Product Code	CSB-PA719421LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCX8
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 21 protein (640-794AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 21 (EC 2.7.11.25) (Mitogen-activated protein kinase kinase kinase MLK4) (Mixed lineage kinase 4), MAP3K21, KIAA1804 MLK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.