



ODR4 Antibody, HRP conjugated

Product Code	CSB-PA725554LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SWX8
Immunogen	Recombinant Human Protein odr-4 homolog protein (153-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein odr-4 homolog (hODR-4) (LAG1-interacting protein) (Transactivated by transforming growth factor beta protein 1), ODR4, C1orf27 TTG1 TTG1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ODR4