



LCE3C Antibody, HRP conjugated

Product Code	CSB-PA725885LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T5A8
Immunogen	Recombinant Human Late cornified envelope protein 3C protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Late cornified envelope protein 3C (Late envelope protein 15) (Small proline-rich-like epidermal differentiation complex protein 3A), LCE3C, LEP15 SPRL3A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LCE3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.