



HES5 Antibody, Biotin conjugated

Product Code	CSB-PA729212LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TA89
Immunogen	Recombinant Human Transcription factor HES-5 protein (1-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor HES-5 (Class B basic helix-loop-helix protein 38) (bHLHb38) (Hairy and enhancer of split 5), HES5, BHLHB38
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HES5