



HES5 Antibody, Biotin conjugated

Product Code	CSB-PA729212LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TA89
Immunogen	Recombinant Human Transcription factor HES-5 protein (1-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor HES-5 (Class B basic helix-loop-helix protein 38) (bHLHb38) (Hairy and enhancer of split 5), HES5, BHLHB38
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HES5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.