



SIRPB1 Antibody, FITC conjugated

Product Code	CSB-PA732939LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TFQ8
Immunogen	Recombinant Human Signal-regulatory protein beta-1 isoform 3 protein (30-206AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal-regulatory protein beta-1 isoform 3 (SIRP-beta-1 isoform 3), SIRPB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SIRPB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.