



MCTP2 Antibody, FITC conjugated

Product Code	CSB-PA738576LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DN12
Immunogen	Recombinant Human Multiple C2 and transmembrane domain-containing protein 2 protein (129-240AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multiple C2 and transmembrane domain-containing protein 2, MCTP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MCTP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.