



ELP2 Antibody, HRP conjugated

Product Code	CSB-PA740224LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IA86
Immunogen	Recombinant Human Elongator complex protein 2 protein (323-457AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 2 (ELP2) (SHINC-2) (STAT3-interacting protein 1) (StIP1), ELP2, STATIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ELP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.