



ELP2 Antibody, HRP conjugated

Product Code	CSB-PA740224LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IA86
Immunogen	Recombinant Human Elongator complex protein 2 protein (323-457AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 2 (ELP2) (SHINC-2) (STAT3-interacting protein 1) (StIP1), ELP2, STATIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ELP2