



IYD Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA754272LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6PHW0 |
| Immunogen | Recombinant Human Iodotyrosine deiodinase 1 protein (24-147AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Iodotyrosine deiodinase 1 (IYD-1) (EC 1.21.1.1) (Iodotyrosine dehalogenase 1), IYD, C6orf71 DEHAL1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | IYD |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |