



IYD Antibody, Biotin conjugated

Product Code	CSB-PA754272LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PHW0
Immunogen	Recombinant Human Iodotyrosine deiodinase 1 protein (24-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iodotyrosine deiodinase 1 (IYD-1) (EC 1.21.1.1) (Iodotyrosine dehalogenase 1), IYD, C6orf71 DEHAL1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IYD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.