



LAMTOR1 Antibody, FITC conjugated

Product Code	CSB-PA757083LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IAA8
Immunogen	Recombinant Human Ragulator complex protein LAMTOR1 protein (2-161AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ragulator complex protein LAMTOR1 (Late endosomal/lysosomal adaptor and MAPK and MTOR activator 1) (Lipid raft adaptor protein p18) (Protein associated with DRMs and endosomes) (p27Kip1-releasing factor from RhoA) (p27RF-Rho), LAMTOR1, C11orf59 PDRO
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LAMTOR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.