



RANBP10 Antibody, HRP conjugated

Product Code	CSB-PA764955LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6VN20
Immunogen	Recombinant Human Ran-binding protein 10 protein (3-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ran-binding protein 10 (RanBP10), RANBP10, KIAA1464
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RANBP10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.