



RUFY3 Antibody, Biotin conjugated

Product Code	CSB-PA765974LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L099
Immunogen	Recombinant Human Protein RUFY3 protein (243-352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein RUFY3 (RUN and FYVE domain-containing protein 3) (Rap2-interacting protein x) (RIPx) (Single axon-regulated protein) (Singar), RUFY3, KIAA0871
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RUFY3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.