



GLRX5 Antibody, HRP conjugated

Product Code	CSB-PA769767LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SX6
Immunogen	Recombinant Human Glutaredoxin-related protein 5, mitochondrial protein (32-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaredoxin-related protein 5, mitochondrial (Monothiol glutaredoxin-5), GLRX5, C14orf87
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLRX5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.