



PTOV1 Antibody, Biotin conjugated

Product Code	CSB-PA769809LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YD1
Immunogen	Recombinant Human Prostate tumor-overexpressed gene 1 protein (2-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostate tumor-overexpressed gene 1 protein (PTOV-1) (Activator interaction domain-containing protein 2), PTOV1, ACID2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PTOV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.