



MEGF8 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA773606LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7Z7M0 |
| Immunogen | Recombinant Human Multiple epidermal growth factor-like domains protein 8 protein (495-663AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Multiple epidermal growth factor-like domains protein 8 (Multiple EGF-like domains protein 8) (Epidermal growth factor-like protein 4) (EGF-like protein 4), MEGF8, C19orf49 EGFL4 KIAA0817 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | MEGF8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |