



MEGF8 Antibody, Biotin conjugated

Product Code	CSB-PA773606LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7M0
Immunogen	Recombinant Human Multiple epidermal growth factor-like domains protein 8 protein (495-663AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multiple epidermal growth factor-like domains protein 8 (Multiple EGF-like domains protein 8) (Epidermal growth factor-like protein 4) (EGF-like protein 4), MEGF8, C19orf49 EGFL4 KIAA0817
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MEGF8