



SH3RF1 Antibody, Biotin conjugated

Product Code	CSB-PA773800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6J0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase SH3RF1 protein (628-742AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase SH3RF1 (EC 2.3.2.27) (Plenty of SH3s) (Protein POSH) (RING finger protein 142) (RING-type E3 ubiquitin transferase SH3RF1) (SH3 domain-containing RING finger protein 1) (SH3 multiple domains protein 2), SH3RF1, KIAA1494 POSH RNF142 SH3MD2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SH3RF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.