



SH3RF1 Antibody, Biotin conjugated

Product Code	CSB-PA773800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6J0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase SH3RF1 protein (628-742AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase SH3RF1 (EC 2.3.2.27) (Plenty of SH3s) (Protein POSH) (RING finger protein 142) (RING-type E3 ubiquitin transferase SH3RF1) (SH3 domain-containing RING finger protein 1) (SH3 multiple domains protein 2), SH3RF1, KIAA1494 POSH RNF142 SH3MD2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SH3RF1