



SFR1 Antibody, Biotin conjugated

Product Code	CSB-PA801256LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86XK3
Immunogen	Recombinant Human Swi5-dependent recombination DNA repair protein 1 homolog protein (5-93AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Swi5-dependent recombination DNA repair protein 1 homolog (Meiosis protein 5 homolog), SFR1, C10orf78 MEI5 MEIR5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SFR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.