



TBC1D1 Antibody, FITC conjugated

Product Code	CSB-PA803114LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TI0
Immunogen	Recombinant Human TBC1 domain family member 1 protein (566-718AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 1, TBC1D1, KIAA1108
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.