



THNSL2 Antibody, HRP conjugated

Product Code	CSB-PA803155LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YJ6
Immunogen	Recombinant Human Threonine synthase-like 2 protein (295-441AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Threonine synthase-like 2 (TSH2) (EC 4.2.3.-) (Secreted osteoclastogenic factor of activated T-cells) (SOFAT), THNSL2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	THNSL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.