



# SLC41A1 Antibody, Biotin conjugated

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA810274LD01HU                                                               |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q8IVJ1                                                                           |
| Immunogen           | Recombinant Human Solute carrier family 41 member 1 protein (1-96AA)             |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Human                                                                            |
| Tested Applications | ELISA                                                                            |
| Form                | Liquid                                                                           |
| Conjugate           | Biotin                                                                           |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified                                                         |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Alias               | Solute carrier family 41 member 1, SLC41A1                                       |
| Immunogen Species   | Homo sapiens (Human)                                                             |
| Research Area       | Signal Transduction                                                              |
| Target Names        | SLC41A1                                                                          |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |