



# ULK2 Antibody, HRP conjugated

Product Code	CSB-PA811629LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYT8
Immunogen	Recombinant Human Serine/threonine-protein kinase ULK2 protein (211-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase ULK2 (EC 2.7.11.1) (Unc-51-like kinase 2), ULK2, KIAA0623
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ULK2